

Metabolomic study of the red algae *Chondrus crispus*

MetaboMER, Metabolomic Platform, FR2424 CNRS-UPMC, Station Biologique de Roscoff, Roscoff, France.

CONTEXT

The goal of this study was to establish the best protocol of extraction of the whole metabolome of algae. Seaweeds contain a lot of polysaccharides. For instance, in red algae, carrageenans confer a gelling nature to tissue, resulting in difficulties to obtain a reproducible manual grinding. Thus, particular attention has to be paid to grinding and automated techniques to carry out such extraction steps, as cryo-grinding, should be preferred. To avoid metabolite degradation, freeze-drying the sample before grinding is important. Thus, we tested automated cryo-grinding for the metabolite extraction of a red algae, *Chondrus crispus*, submitted to elicitation by methyl jasmonate (MeJA).

MATERIAL

- Precellys 24 homogenizer with Cryolys cooling device to have a constant temperature of -20 °C within the homogenization chamber using liquid nitrogen.
- Precellys lysing kit : 03961-1-008 (MK28-R).
- Sample: 3x3 replicates of MeJA elicited & control *Chondrus crispus*.
- Extraction solvent: 1 ml ice-cold [MeOH/H₂O (8:2)] for 200mg fresh tissue (added after homogenization).

PROTOCOL

- Each lysing tube was added of 150 mg Fontainebleau sand.
- Precellys 24: 1x30s 6800rpm - 1x50s 6000rpm - 1x30s 6800rpm - 1x90s 6200rpm (with 5min break between each run).
- Extraction was carried out at +4°C during 1hour with constant agitation. Centrifugation 3000g, 15 min.
- Samples were injected in UHPLC/MS on Dionex U3000 RSLC coupled to a LTQ-Orbitrap Discovery (Thermo, Les Ulis). Data treatment with XCMS and Mev.

RESULTS

Our results clearly demonstrate the elicitation of the red algae *Chondrus crispus* by MeJA (results not shown). Moreover, the 9 different replicates display very reproducible metabolomic profiles (Figure 1), thanks to simultaneous grinding with Precellys 24. During this study, the Cryolys guaranteed a constant temperature of -20°C within the homogenization chamber and thus avoided an uncontrolled defrosting of the frozen samples.





Figure 1: Hierarchical clustering of all detected ions from control and MeJA elicited *Chondrus crispus* samples. Samples are in columns; each line represents a detected ion and its variations, from low amount (green) to higher amount (red). Here, only the 91 first ions are represented (2448 detected).



CNRS UPMC

Station Biologique
Roscoff



03712-810-DU068



CONCLUSION

Compare to manual grinding, simultaneous cryo-grinding avoided biological material degradation and allowed higher reproducibility, faster grinding and better metabolite extraction. Cryo-grinding with the combo Precellys 24 and Cryolys is for our platform a daily necessary tool to perform robust and high-throughput metabolomic analysis.

For more details, please contact precellys@bertin.fr

