

# Xanthones purification with Centrifugal Partition Chromatography SCPC-250

# Application Note 02

### Introduction

The purple mangosteen (garcinia mangostana) is a tree of clusiaceae family. Fruits and leaves are used in medicine. **Centrifugal Partition Chromatography (CPC)** also known as **Counter Current Chromatography (CCC)** is a preparative, pilot and industrial liquid purification technique that does not require traditional solid supports. CPC was used to purify few mg of xanthones from an extract of garcinia mangostana pericarp.



# Materials and Methods

A Gilson <u>SCPC-250</u> coupled with a Gilson <u>PLC2250</u> system and Lachrome Elite HPLC with DAD and ESI-MS detectors were used.

CPC solvent system is determined with shake flask method to get a Kd= [HPLC peak area of anthocyanin] <sub>stat</sub>/ [HPLC peak area of anthocyanin] <sub>mobile</sub> closed to one.







### Table 1: CPC conditions

CPC column volume	: 250 mL
Flow rate	: 4 mL/min
Rotation speed	: 1700 rpm
Solvent system	: Arizona R
Mass injected	: 175 mg of ethanolic extract
Detection	: 320 and 254 nm

### **Results and Discussion**

During CPC purification, automated HPLC injections are performed (in red on figure 1).

Three major compounds of the extract were isolated and there structures were identified by LC-MS :

- $3-isomangostin [(M+H+H_2O)]^+ = 429]$
- gartanin [(M+H)<sup>+</sup> = 397]
- $\alpha$ -mangostin [(M+H)<sup>+</sup> = 411]



Figure 1: CPC chromatogram at 254 nm (in black) and 320 nm (in orange)

Gilson www.gilson.com/asia asia@gilson.com





# Conclusions

The combination of CPC, PLC and HPLC-ESI-MS allows the fast and simultaneous separation and identification of natural xanthones in a crude extract whatever the mobile phase used for the separation.

Interfacing PLC/CPC with mass spectrometry provides a new automated analytical methodology in screening crude natural extracts. It's an isolation procedure that combines the advantages of CPC with the low detection limit of mass spectrometry.

