



Analysis of MAb and MAb Fragments on Antibodix™ WCX NP5

Ion Exchange Chromatography is frequently used for antibody analysis. Antibodies and antibody fragments can all be separated on weak cation exchange columns based on their charge states. Antibodix™ WCX columns are specially designed for high resolution, high efficiency and high recovery separations of antibodies.


In this study, we investigated antibody fragments such as heavy and light chains, Fab/Fc and F(ab')₂ using Antibodix™ WCX separation.

Highlighted FACTS:


- ▶ Antibodix™ WCX NP5 4.6 x 250 mm can successfully separate monoclonal antibody variants under different mobile phase systems such as pH and salt gradients.
- ▶ Monoclonal antibody purity, heterogeneity and stability can be monitored using Antibodix™ WCX NP5.
- ▶ The smaller particle size in Antibodix™ WCX NP5 offers superior resolution in comparison to Antibodix™ WCX NP10.

MAb loading test on Antibodix™ WCX NP5

Column: Antibodix™ WCX NP5 (4.6x250 mm)
 Mobile phase: A: 20 mM Phosphate buffer pH 5; B: A+10 mM NaCl pH 7.5
 Flow rate: 0.8 mL/min; Temperature: 30 °C;
 Sample: Intact MAb 321 (5 mg/mL)




MAb analysis on Antibodix™ NP5 and NP10



Column: Antibodix™ WCX NP5 and Antibodix™ WCX NP10 (4.6x250 mm)
 Mobile phase: A: 20mM Sodium Acetate pH 5.15, B: A + 1M LiCl, Flow rate: 0.8 mL/min;
 Injection volume: 20 µL,
 Sample: MAb 321 (5 mg/mL)


Analysis of MAB 321 using a LiCl gradient

Column: Antibodix™ WCX NP5 (4.6x250 mm)
 Mobile phase: A: 20 mM Sodium Acetate pH 5.15, B: A + 1 M LiCl
 Flow rate: 0.8 mL/min; Injection volume: 20 µL, Detection: UV280 nm,
 Temperature: 30 °C; Sample: Intact MAb 321 (5 mg/mL)




Analysis of MAB 321 using a NaCl and pH gradient

Column: Antibodix™ WCX NP5 (4.6x250 mm)
 Mobile phase: A: 20 mM Phosphate buffer pH 5; B: A+10 mM NaCl pH 7.5
 Flow rate: 0.8 mL/min; Injection volume: 20 µL; Detection: UV280 nm
 Temperature: 30 °C; Sample: Intact MAB 321 (5 mg/mL)



MAb stability test on Antibodix™ WCX NP5

Column: Antibodix™ WCX NP5
 Mobile phase: A: 20 mM Sodium Acetate pH 5.15, B: A + 1 M LiCl, Flow rate: 0.8 mL/min;
 Injection volume: 20 µL,
 Sample: Intact MAB 321 (5 mg/mL)





Analysis of MAb and MAb Fragments on Antibodix™ WCX NP5

What is Antibodix™ WCX NP5


Antibodix™ WCX NP5 (Weak Cation Exchange):

Comprised of rigid, spherical, highly cross-linked non-porous poly(styrene divinylbenzene) (PS/DVB) beads. The PS/DVB particle surface is grafted with a hydrophilic, neutral polymer layer which is nanometers thick. The resin surface is covered by a hydrophilic coating which eliminates non-specific bindings with antibody proteins, leading to high efficiency and high recovery separations. On top of the hydrophilic layer, weak cation-exchange functional groups are attached via a proprietary chemistry, resulting in a high capacity ion-exchange layer.

Technical Specifications:

Phase	Antibodix™ WCX NP5
Material	Weak cation exchange groups bonded to a hydrophilic film grafted on PS/DVB
Particle size (µm)	5
Pore size (Å)	Non-porous
pH stability	2 – 12
Backpressure (psi)	~ 2,500
Maximum backpressure	~ 6,000
Maximum temperature	~ 80 °C
Mobile phase compatibility	Aqueous or a mixture of water and acetonitrile, acetone or methanol


Comparison of Sepax's Antibodix™ NP5 to a competitor's column



Column: Antibodix™ WCX NP5 (4.6x250 mm)
 Mobile phase:
 A: 20 mM Sodium Acetate pH 5.15,
 B: A + 1 M LiCl
 Flow rate: 0.8 mL/min;
 Temperature: 30 °C;
 Sample: Intact MAb 321 (5 mg/mL)
 Injection volume: 20 µL


Analysis of Antibody Fragments on Antibodix™ WCX NP5

Papain digested MAb analysis on Antibodix™ NP5 and SDS page gel




Column: Antibodix™ WCX NP5 (4.6x250 mm)
 Mobile phase:
 A: 20mM acetic acid+50mM NaCl pH3.5
 B: 20mM sodium succinate+50mM NaCl pH6
 Sample: 100µg Papain digest MAb321
 Gradient: 5 min 30% B, 0.8 mL/min;
 25 min 85%-100% B, 0.65mL/min

Fab/Fc loading test on Antibodix™ NP5 4.6 x 250 mm




Column: Antibodix™ WCX NP5 (4.6x250 mm)
 Mobile phase:
 A: 20mM acetic acid+50mM NaCl pH3.5
 B: 20mM sodium succinate+50mM NaCl pH6
 Sample: Papain digest MAb321
 Gradient: 5 min 30% B, 0.8 mL/min;
 25 min 85%-100% B, 0.65mL/min

Fab/Fc on Antibodix™ NP5



Column: Antibodix™ WCX NP5 (4.6x250 mm)
 Mobile phase:
 A: 20 mM Phosphate Buffer pH 5.5
 B: A + 1 M NaCl
 Flow rate: 0.8 mL/min;
 Detection: UV 280 nm
 Sample: 25 µg Papain digested MAB

F(ab')₂ on Antibodix™ NP5



Column: Antibodix™ WCX NP5 (4.6x250 mm)
 Mobile phase:
 A: 20 mM phosphate buffer, pH 5.5,
 B: A+ 1 M NaCl
 Flow rate: 0.8 mL/min;
 Detection: UV 280 nm
 Sample: 50 µg Pepsin digested MAB

Sepax Technologies, Inc.

5 Innovation Way

Newark, Delaware 19711, USA

Tel: (302) 366-1101 | Fax: (302) 366-1151

E-mail: info@sepax-tech.com