



RNA EXTRACTION FROM NORMAL DOG SKIN BIOPSIES

Dr Nicky Craig

College of Veterinary Medicine, University of Florida, Gainesville, Florida, USA

/ CONTEXT

Skin is very fibrous and elastic, and efficient RNA extraction from this tissue can be very difficult. We have previously tested several homogenizers to homogenize dog skin samples, with unsatisfactory results, including poor homogenization, poor RNA yield and warming of the samples. To find a better option for RNA extraction from dog skin samples, several types of beads were tested for efficacy in homogenizing frozen skin biopsies using a Precellys® 24 bead homogenizer.

Lysis buffer in the tubes was defrosted on ice and the samples homogenized by running 2-3 cycles of 15 sec at 5,000 rpm using a Precellys® 24. Whole frozen biopsies in 500µl of ice-cold lysis buffer were also homogenized for comparison.

RNA Extraction: RNA was purified using a Qiagen RNeasy kit (Trizol and RLT samples) or 5Prime PerfectPure RNA Cell kit, in a final volume of 60µl. RNA quality was analysed using Nanodrop 2000 +/- Agilent TapeStation analysis (RNA integrity number, RIN).

/ MATERIALS

Tissue Samples: Adjacent 4mm full thickness skin biopsies taken from normal 9 week old dog cadaver skin which had been frozen immediately post-mortem and stored at -80°C.

Homogenizer: Precellys® 24.

Lysing Beads tested:

- 1mm zirconium oxide beads (Benchmark Scientific in MP Biomedicals 2ml tubes).
- 2mm zirconium oxide beads (Benchmark Scientific in MP Biomedicals 2ml tubes).
- 2.8mm and 5.0mm zirconium oxide mix beads (CKMix50_2mL, Bertin, Cat #: KT03961-1-013.2).
- 0.7mm garnet flakes and one 6mm zirconium oxide bead (GK60_2mL, Bertin, Cat #: KT03961-1-015.2).
- 6.35mm zirconium oxide beads (Lysing Matrix M, MP Biomedicals).
- One or two 6.35mm steel beads (Lysing Matrix SS, MP Biomedicals).

Lysis Buffers tested: Trizol, RLT+beta mercaptoethanol (Qiagen) or PerfectPure RNA Cell kit Lysis Buffer+TCEP (5Prime).

/ RESULTS

Dicing the tissue prior to homogenization was essential as whole 4mm biopsies didn't homogenize well.

There did not appear to be a linear relationship between RNA concentration and bead size or density. It is likely the variation in RNA concentration observed is mostly due to variation in how finely diced the tissue was before homogenization.

Trizol has previously been observed to turn black in combination with steel beads. This didn't occur here, possibly because the sample didn't get warm during homogenization.

Table1. RNA concentrations and purity observed using different lysing kits.

Homogenization (diced tissues)	RNA ng/µl	A 260/280	A 260/230	RIN
Trizol ZrO ₂ 1mm, BS in MP Bio tubes	293.4	2.04	2.17	5.5
Trizol ZrO ₂ 2mm, BS in MP Bio tubes	247.3	2	2.07	
Trizol garnet flakes, ZrO ₂ 6mm, Bertin	110.8	2.01	1.9	
Trizol ZrO ₂ 6.35mm, MP Bio	263.2	2.03	2.18	
Trizol steel 6.35mm x1, MP Bio	372.9	2.04	2.18	
Trizol steel 6.35mm x2, MP Bio	346.4	2.04	2.19	5.5
RLT ZrO ₂ 1mm, BS in MP Bio tubes	284.8	2.04	2.14	7.1
RLT ZrO ₂ 2mm, BS in MP Bio tubes	218.4	2.03	2.04	7.8
RLT ZrO₂ 2.8mm & 5mm, Bertin	371.1	2.04	2.16	8.2
RLT garnet flakes, ZrO ₂ 6mm, Bertin	222.9	2.04	2.21	
RLT ZrO ₂ 6.35mm, MP Bio	198	2.04	2.16	
RLT steel 6.35mm x1, MP Bio	269	2.04	2.16	8.6
PerfectPure ZrO ₂ 1mm, BS in MP Bio tubes	203	2.09	2.2	7.3
PerfectPure ZrO₂ 2.8mm & 5mm, Bertin	369	2.07	2.15	8.3
PerfectPure steel 6.35mm, MP Bio	349	2.07	1.86	4.1

/ PROTOCOL

Sample preparation: Beads in 2ml tubes were rinsed with 200µl of lysis buffer, then 400µl of lysis buffer was added to the tube and it was placed on dry ice. A 50µl droplet of lysis buffer was frozen in a plastic petri dish placed on a block of dry ice. Frozen 4mm biopsies were placed on the frozen lysis buffer and another 50µl of ice-cold lysis buffer was slowly placed on top and allowed to freeze.

A sterile scalpel blade was used to dice the frozen tissue and lysis buffer. Diced tissue and lysis buffer were added to the frozen beads and lysis buffer in the 2ml tubes, and the tubes were stored at -80°C until homogenization.

/ CONCLUSION

The Precellys® 24 tissue homogenizer efficiently homogenized finely diced dog skin within a short time (30-45 seconds) without warming the samples. Higher concentrations and better quality RNA were purified using the Precellys® 24 when compared to the bead and rotor-stator homogenizers we have tested previously. Some samples exhibited a degree of RNA degradation, as indicated by a RIN of less than 8 on TapeStation analysis.

SPPP-810-DU131



□東京 〒162-0805 東京都新宿区矢来町 113 番地 TEL (03)3235-0661(代) / FAX (03)3235-0669
 □大阪 〒532-0005 大阪市淀川区三国本町2丁目12番4号 TEL (06)6396-0501(代) / FAX (06)6396-0508
 □福岡 〒812-0054 福岡市東区馬出 1 丁目 2 番 23 号 TEL (092)631-1012(代) / FAX (092)641-1285

※会社名および商品名は、各会社の商標または登録商標です。
 ※本カタログに記載の規格・仕様・外観は予告なく変更する場合がありますので御了承下さい。

<http://www.technosaurus.co.jp>

【Precellys を用いたイヌ皮膚からの RNA 抽出】 Precellys-APS-1710-93