

Maintaining yeast viability during tissue homogenization with Minilys

Lab. Microbiologie Fondamentale et Pathogénicité CNRS UMR-5234 - Unité Candida et Pathogénicité, Bordeaux, France

CONTEXT

Fungi of medical interest in recent years have become primarily responsible for opportunist parasitic diseases, mainly because of the strong increase of immunodepressed patients. Fungal infections pose a challenge to health professionals who have to face variety of new emerging pathogens, but also manage the problem of resistance to antifungal treatments. The team focus its research on yeasts of medical interest belonging to the genus *Candida*.

MATERIAL

- Minilys homogenizer.
- Precellys kits: CK28R_2mL (03961-1-007) + 2CK50 (03961-1-106); CK28_7mL (03961-1-302) + 4CK50 (03961-1-106);
- Samples: Duplicate of whole tissues of mice (Lung, heart, liver, brain, kidney and spleen).
- Buffer: 1mL water in 2 mL vial or 2mL water in 7mL.
- *Candida* culture: 1.10⁷ CFU; 100µL in 2mL vial or 200µL in 7 mL vial.

PROTOCOL

- Minilys setting: 5000 rpm.

	Lung	Heart	Liver	Brain	Kidney	Spleen
	[0,63g / 0,66g]	[0,110g / 0,02g]	[1,18g / 1,10g]	[0,35g / 0,41g]	[0,43g / 0,45g]	[0,33g / 0,22g]
2mL	2x60s	2x30s				
7mL			2x15s	10s	30s	30s

- Analysis: Assessment of *Candida* concentration (CFU/mL) by cultural method after spreading 100µL of a 100 dilution in YPD plates.

RESULTS

Assessment of yeast concentration was performed on the homogenate samples (Results not shown).

Effect of homogenization on the viability of yeast cells was performed on a pure *Candida* culture at 1.10⁷ CFU. Following Minilys homogenization, either performed in 7 mL or 2mL vial, 100 % of yeast cells survived after the maximum of time set independently the tissue. It's important to limit the temperature rise during long homogenization (>40sec.) by keeping the samples in ice during a few minutes between 40s cycles (Fig.1).

Effect on viability of yeast cells




Figure 1: Effect of tissue homogenization on the viability of *Candida*. Data represents the % of yeast survived (CFU/ml) after the duration of homogenization specified.

CONCLUSION

Homogenization of whole tissues of mice (Lung, heart, liver, brain, kidney and spleen) using the **Minilys** with **Precellys lysing kits** does not affect the viability of the yeast pathogen *Candida*. The yeast viability is maintained during tissue homogenization up to 2 min limiting excessive temperature rise.

The potential to process individual samples in a short time period with no risk of cross-contamination has considerable benefits over traditional homogenization methods.

For more details, please contact precellys@bertin.fr

03712-810-DU072



イムス機器株式会社

<http://www.technosaurus.co.jp>

- 東京
〒162-0805 東京都新宿区矢来町 113 番地 TEL(03)3235-0661(代) FAX(03)3235-0669
- 大阪
〒532-0005 大阪市淀川区三国本町2丁目12番4号 TEL(06)6396-0501(代) FAX(06)6395-2588
- 福岡
〒812-0054 福岡市東区馬出 1 丁目 2 番 23 号 TEL(092)631-1012(代) FAX(092)641-1285

※会社名および商品名は、各会社の商標または登録商標です。
 ※本カタログに記載の規格・仕様・外観は予告なく変更する場合がありますので御了承下さい。