

# Proteomix™ SCX 1.7 μm - 4.6 x 50 mm

Column: Proteomix™ SCX-NP1.7 ( 1.7 μm, 4.6 x 50 mm)

Mobile phase: A: 10 mM Phosphate buffer (pH 7.0), B: A + 0.5 M NaCl

Flow rate: 0.6 mL/min, Detector: UV 280 nm, Column temperature: 25°C

Injection volume: 5 μL, Samples: 5 mg/mL Erbitux and biosimilar candidate, Pressure: 230 bar

